

LIEBHERR LOADING PLAN

SCIENTIFIC REFRIGERATOR LKv 5710




LOADING WITH PFIZER/BIONTECH TRAYS FOR 195 VIALS

General Information

The temperature-controlled Pfizer Cool Box enables the logistic of the Pfizer / BioNTech vaccine BNT162b2 at a temperature of -70°C by usage of dry-ice. One box/container can be filled with 5 trays, each tray with 195 vials. According to available information 5 vaccination doses can be taken from one vial.

The reusable Pfizer Cool Box consists of an outer carton and a payload sleeve for the five trays. The gap between outer carton and payload sleeve is filled up with dry-ice. With one dry-ice-loading the container keeps a temperature of -70°C for at least 10 days. The vaccines can be stored at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$ for max. 5 days before used.



Vaccine boxes / Payload (5 trays): height 7.9" = 20.1 cm x width 9.1" = 23.2 cm x depth 9.1" = 23.2 cm


Transport container / outer carton: height 22" = 55.9 cm x width 15.7" = 39.9 cm x depth 15.7" = 39.9 cm

Source:
https://si.wsj.net/public/resources/images/E3-HG530_COVIDVA_300PX_202009_03213118.jpg


LIEBHERR

LKv 5710 Loading Capacity Pfizer / BioNTech BNT162b2

24 vaccine-boxes, each vaccine box with 5 trays + 2 additional trays (green marked in sketch), **overall 122 trays**;
Each tray with 195 vials for 5 doses = **23.790 vials / 118.950 doses** (or people that can get a vaccination)



LKv 5710



Disclaimer: All information regarding the Pfizer Cool Box as well as its dimensions are without guarantee. This information refers to existing press reports and may change due to further developments.